


Introduction to FreeSurfer

<http://surfer.nmr.mgh.harvard.edu>

Allison Stevens

Bruce Fischl, Doug Greve, Nick Schmansky, Jenni Pacheco
freesurfer@nmr.mgh.harvard.edu



Cortical Surface Reconstruction

FreeSurfer creates computerized models of the brain from MRI data.

Input:

- T1-weighted (MPRAGE,SPGR)
- 1mm³ resolution


Cortical Surface Reconstruction

- Finds white/gray boundary – wm surface
- Finds pial/CSF boundary – pial surface
- To “Find” uses:
 - Intensity information, spatial location, geometric structure
 - Tessellation, neighbors, talairach coordinates
- Subcortical Segmentation

Surface Model



- Mesh (“Finite Element”)
- Vertex = point of 6 triangles
- XYZ at each vertex
- Triangles/Faces ~ 150,000
- Area, Distance
- Curvature, Thickness



Cortical Reconstruction Goals

- Geometrically Accurate surfaces
 - Accurately follow the boundaries seen on the scan for each of your individual subjects
- Topologically Correct surfaces
 - Each surface is a 2-D continuous, non self-intersecting sheet and can be inflated into a perfect sphere
- Surfaces are only as good as your scan.


MR Anatomy Caveats

- Dependent on data quality
 - Contrast to noise
 - Signal to noise
 - Voxel resolution
- MR Artifacts
 - MR susceptibility
 - MR distortions
- Variations in MR tissue parameters across regions of the brain are altered in different populations


FreeSurfer Output

- Volumes
- Surfaces
- Surface Overlays
- ROI Summaries


Volumes



orig.mgz



T1.mgz



brainmask.mgz



wm.mgz




filled.mgz

Subcortical Mass



aseg.mgz




aparc+aseg.mgz


- `$$SUBJECTS_DIR/bert/mri`
- All “Conformed” 256^3 , 1mm

Volume Viewer:
tkmedit


Surfaces




orig




white




pial



inflated



sphere,sphere.reg




flat

- `$SUBJECTS_DIR/bert/surf`
- Number/Identity of vertices stays the same (except flat)
- XYZ Location Changes
- Flattening not done as part of standard reconstruction


Surface Viewer:
tksurfer

Surface Overlays


lh.sulc on inflated



lh.curv on inflated



lh.thickness on inflated





lh.sulc on pial



lh.curv on inflated



fMRI on flat



lh.aparc.annot on inflated

- Value for each vertex
- Color indicates value
- Color: gray, red/green, heat, color table
- Rendered on any surface
- fMRI/Stat Maps too

ROI Summaries

aseg.stats

- volumes of subcortical structures (mm^3)

aparc.stats

- thickness of cortical parcellation structures (mm)
- total white matter volume (mm^3)
- number of vertices in cortex
- surface area of cortex (mm^2)

ROI Summaries: Make Your Own

yourROI.label

- Draw your own surface label
- use `mris_anatomical_stats` to get data

mris_volume

- Total volume within a surface you specify

mris_wm_volume

- Total volume within white surface ignoring non-wm voxels in `aseg.mgz`


Reconstruction Environment

- Installation directory:
 \$FREESURFER_HOME
- Set-up Environmental Variables
- Unix command-line (Linux, MacOSX)
- Directory structure, naming conventions
- File Formats

Set-up Environmental Variables

Subject ID


`$SUBJECTS_DIR`



Before running FreeSurfer, **must** set `$FREESURFER_HOME` and `$SUBJECTS_DIR`


FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)




FreeSurfer Directory Tree

Directories used often are in green.



MGZ File Format



- mgz = compressed MGH file
- Can store 4D (like NIFTI)
- cols, rows, slices, frames
- Generic: volumes and surfaces

- Eg, Typical Anatomical volume: 256 x 256 x 128 x 1



- FreeSurfer can read/write:
NIFTI, Analyze, MINC
Careful with NIFTI! (32k
limit)

- FreeSurfer can read:
DICOM, Siemens
IMA, AFNI

Other FreeSurfer File Formats

Unique to FreeSurfer

- Surface: lh.white, lh.pial, lh.orig
- Curv: lh.curv, lh.sulc, lh.thickness
- Annotation: lh.aparc.annot
- Label: lh.pericalcarine.label

Starting the Reconstruction Process


Before running FreeSurfer, **must** set
\$FREESURFER_HOME and \$SUBJECTS_DIR

**recon-all -i /path/to/your/raw/data1 -i
/path/to/your/raw/data2 -all -s subject_id**

- This will create the subject directory ‘**subject_id**’ in your \$SUBJECTS_DIR and convert your **2 raw acquisitions** to mgz and use them as input for the ‘**-all**’ command.

Alternative: Add Your Data

- `cd $SUBJECTS_DIR`
- `mkdir -p bert/mri/orig`
- `mri_convert yourdicom.dcm bert/mri/orig/001.mgz`
- `mri_convert yourdicom.dcm bert/mri/orig/002.mgz`
- `recon-all -all -s bert`



Individual Steps

Volumetric Processing Stages (subjid/mri):

1. Motion Cor, Avg, Conform (orig.mgz)
2. Non-uniform inorm (nu.mgz)
3. Talairach transform computation
(talairach/talairach.xfm)
4. Intensity Normalization 1 (T1.mgz)
5. Skull Strip (brainmask.mgz)
6. EM Register (linear volumetric registration)
7. CA Intensity Normalization (norm.mgz)
8. CA Non-linear Volumetric Registration
9. CA Label (Volumetric Labeling) (aseg.mgz)
10. Intensity Normalization 2 (T1.mgz)
11. White matter segmentation (wm.mgz)
12. Edit WM With ASeg
13. Fill and cut (filled.mgz)

Blue = Manual Intervention

recon-all -help

Surface Processing Stages (subjid/surf):

14. Tessellate (?h.orig.nofix)
15. Smooth1
16. Inflate1
17. QSphere (?h.qsphere)
18. Automatic Topology Fixer (?h.orig)
19. Final Surfs (?h.white ?h.pial ?h.thickness)
20. Smooth2 (?h.smoothwm)
21. Inflate2 (?h.inflated)
22. Aseg Statistics (stats/aseg.stats)
23. Cortical Ribbon Mask (?h.ribbon.mgz)
24. Spherical Morph
25. Spherical Registration (?h.sphere.reg)
26. Map average curvature to subject
27. Cortical Parcellation (Labeling)
28. Cortical Parcellation Statistics
29. Cortical Parcellation mapped to Aseg
30. White Matter Parcellation (wmparc.mgz)

Note: ?h.orig means lh.orig or rh.orig

Reconstruction Stages

recon-all is broken into three stages

- autorecon1
- autorecon2
- autorecon3

-autorecon1

Volumetric Processing Stages (subjid/mri):


1. Motion Cor, Avg, Conform (orig.mgz)
2. Non-uniform inorm (nu.mgz)
3. Talairach transform computation
(talairach/talairach.xfm)
4. Intensity Normalization 1 (T1.mgz)
5. Skull Strip (brainmask.mgz)
6. EM Register (linear volumetric registration)
7. CA Intensity Normalization (norm.mgz)
8. CA Non-linear Volumetric Registration
9. CA Label (Volumetric Labeling) (aseg.mgz)
10. Intensity Normalization 2 (T1.mgz)
11. White matter segmentation (wm.mgz)
12. Edit WM With ASeg
13. Fill and cut (filled.mgz)

Surface Processing Stages (subjid/surf):


14. Tessellate (?h.orig.nofix)
15. Smooth1
16. Inflate1
17. QSphere (?h.qsphere)
18. Automatic Topology Fixer (?h.orig)
19. Final Surfs (?h.white ?h.pial ?h.thickness)
20. Smooth2 (?h.smoothwm)
21. Inflate2 (?h.inflated)
22. Aseg Statistics (stats/aseg.stats)
23. Cortical Ribbon Mask (?h.ribbon.mgz)
24. Spherical Morph
25. Spherical Registration (?h.sphere.reg)
26. Map average curvature to subject
27. Cortical Parcellation (Labeling)
28. Cortical Parcellation Statistics
29. Cortical Parcellation mapped to Aseg
30. White Matter Parcellation (wmparc.mgz)

recon-all -help


Motion Correction and Averaging



Does not change native resolution.



Conform




Changes to 256^3 , 1mm^3
All volumes will be conformed.



orig Volume

Non-Uniform Intensity Correction


- Uses MNI tool
- Removes B1 bias field



nu Volume

Talairach Transform

- Computes 12 DOF transform matrix
- Does NOT resample
- MNI305 template
- Used to help find structures (eg, CC)
- Can also be used to localize functional activation
- `mri/transforms/talairach.xfm`



Intensity Normalization


- Presegmentation (T1.mgz)
 - All WM = 110 intensity
 - Pre- and Post-Skull Strip



T1 Volume


Skull Strip

- Removes all non-brain
 - Skull, Eyes, Neck, Dura
- brainmask.mgz



FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)



-autorecon2

Volumetric Processing Stages (subjid/mri):

1. Motion Cor, Avg, Conform (orig.mgz)
2. Non-uniform inorm (nu.mgz)
3. Talairach transform computation
(talairach/talairach.xfm)
4. Intensity Normalization 1 (T1.mgz)
5. Skull Strip (brainmask.mgz)
6. EM Register (linear volumetric registration)
7. CA Intensity Normalization (norm.mgz)
8. CA Non-linear Volumetric Registration
9. CA Label (Volumetric Labeling) (aseg.mgz)
10. Intensity Normalization 2 (T1.mgz)
11. White matter segmentation (wm.mgz)
12. Edit WM With ASeg
13. Fill and cut (filled.mgz)


Surface Processing Stages (subjid/surf):

14. Tessellate (?h.orig.nofix)
15. Smooth1
16. Inflate1
17. QSphere (?h.qsphere)
18. Automatic Topology Fixer (?h.orig)
19. Final Surfs (?h.white ?h.pial ?.thickness)
20. Smooth2 (?h.smoothwm)
21. Inflate2 (?h.inflated)
22. Aseg Statistics (stats/aseg.stats)
23. Cortical Ribbon Mask (?h.ribbon.mgz)
24. Spherical Morph
25. Spherical Registration (?h.sphere.reg)
26. Map average curvature to subject
27. Cortical Parcellation (Labeling)
28. Cortical Parcellation Statistics
29. Cortical Parcellation mapped to Aseg
30. White Matter Parcellation (wmparc.mgz)

recon-all -help

Note: lh processed completely first, then rh.



Automatic Volume Labeling



ASeg Volume


- Used to determine volumes of subcortical structures
- Used to fill in subcortical structures for creating subcortical mass
- aseg.mgz

Volume-based Labeling



Labeling is determined by location and intensity.


Validation of Volume Labeling *



Manual labeling done by Center for Morphometric Analysis (CMA)

*Thanks to Drs Larry Seidman and Jill Goldstein for providing this data.


Find “Subcortical Mass”



- All White Matter
- All Subcortical Structures
- Ventricles
- Excludes brain stem and cerebellum
- Hemispheres separated
- Completely connected (no islands)
- Many Stages ... More Later ...

White Matter Segmentation

- Separates white matter from everything else
- “Fills in” subcortical structures
- Cerebellum removed, brain stem still there




mri_segment

mri_edit_wm_with_aseg

mri_preless


Fill and Cut (Subcortical Mass)

- Fills in any voids
- Removes any islands
- Removes brain stem
- Separates hemispheres (each hemi has different value)
- filled.mgz = “Subcortical Mass”





FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)




Tessellation




orig surface
surf/lh.orig
surf/rh.orig

- Mosaic of triangles (“tessellation”)
- Errors: Donut holes, handles
 - Subsequently fixed by the automatic topology fixer




Inflation: Visualization




Dale and Sereno, 1993; Dale et al., 1999; Fischl et al., 1999; Fischl et al., 2000; Fischl et al., 2001


Automatic Topology Fixer




White Matter Surface




- Nudge orig surface
- Follow T1 intensity gradients
- Smoothness constraint
- Vertex Identity stays constant




Pial Surface




- Nudge white surface
- Follow T1 intensity gradients
- Vertex Identity Stays




Optimal Surface Placement




Gray-White Boundary




Outer Cortical Surface




Gray/CSF Deformation



Gray-White Boundary




Outer Cortical Surface



Dale and Sereno, 1993; Dale et al., Dale et al., 1999; Fischl et al., 1999;
Fischl et al., 2000; Fischl et al., 2001


Cortical Thickness




Thickness Maps

- Red regions are thinner
- Yellow regions are thicker


18M




48M




88M




18F



44F



88F







1mm

2mm

3mm




Histological Validation



FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)



-autorecon3

Volumetric Processing Stages (subjid/mri):

1. Motion Cor, Avg, Conform (orig.mgz)
2. Non-uniform inorm (nu.mgz)
3. Talairach transform computation
(talairach/talairach.xfm)
4. Intensity Normalization 1 (T1.mgz)
5. Skull Strip (brainmask.mgz)
6. EM Register (linear volumetric registration)
7. CA Intensity Normalization (norm.mgz)
8. CA Non-linear Volumetric Registration
9. CA Label (Volumetric Labeling) (aseg.mgz)
10. Intensity Normalization 2 (T1.mgz)
11. White matter segmentation (wm.mgz)
12. Edit WM With ASeg
13. Fill and cut (filled.mgz)


Surface Processing Stages (subjid/surf):

14. Tessellate (?h.orig.nofix)
15. Smooth1
16. Inflate1
17. QSphere (?h.qsphere)
18. Automatic Topology Fixer (?h.orig)
19. Final Surfs (?h.white ?h.pial ?h.thickness)
20. Smooth2 (?h.smoothwm)
21. Inflate2 (?h.inflated)
22. Aseg Statistics (stats/aseg.stats)
23. Cortical Ribbon Mask (?h.ribbon.mgz)
24. Spherical Morph
25. Spherical Registration (?h.sphere.reg)
26. Map average curvature to subject
27. Cortical Parcellation (Labeling)
28. Cortical Parcellation Statistics
29. Cortical Parcellation mapped to Aseg
30. White Matter Parcellation (wmparc.mgz)


recon-all -help

Note: lh processed completely first, then rh.

Surface-Based Spherical Coord System





“Spherical” Registration



Sulcal Map

Spherical Inflation




High-Dimensional
Registration to
Spherical Template




Template: average.curvature.filled.buckner.40.tiff


Spherical Inflation




Inflated Surface




Transformed Surface




Spherical Registration to Atlas




Individual Subject



Atlas (Target)





Cortical Parcellation



Spherical Template
based on Manual
Parcellation




Map to Individual
Thru Spherical Reg



Fine-tune based on
individual anatomy

Cortical Parcellation: 2




Atlases: [curvature.buckner40filled.desikan_killiany](#), [atlas_2005_simple](#)

Thanks to Christophe Destrieux for this slide.


FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)



FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)



Actual Workflow: version 1

1. recon-all -all (Stages 1-20) ~30 hours
2. Check talairach transform, skull strip, normalization
3. Check surfaces
 1. Add control points: recon-all -autorecon2-cp -autorecon3 (Stages 10-29)
 2. Edit wm.mgz: recon-all -autorecon2-wm -autorecon3 (Stages 13-29)
 3. Edit brain.mgz: recon-all -autorecon2-pial (Stage 19-22)

Note: all stages can be run individually

Actual Workflow: version 2

1. recon-all --autorecon1 (Stages 1-5) ~45 min
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22) ~20 hours
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29) ~6 hours

Note: all stages can be run individually

Tutorials

<https://surfer.nmr.mgh.harvard.edu/fswiki/FsTutorial>

On Linux Machines:

- `setenv` `FREESURFER_HOME` `/usr/pubsw/freesurfer`
- `source` `$FREESURFER_HOME/SetUpFreeSurfer.csh`
- `setenv` `SUBJECTS_DIR` `/ircuser/FSWorkshop/buckner_data/tutorial_subjs`

On Macs:


- `setenv` `FREESURFER_HOME` `/Applications/freesurfer`
- `source` `$FREESURFER_HOME/SetUpFreeSurfer.csh`
- `setenv` `SUBJECTS_DIR` `/FSWorkshop/buckner_data/tutorial_subjs`
- Use the volume and surface viewing tools to look at correct output from all steps.

Troubleshooting

- Skull Strip Errors
- Intensity Normalization
- Segmentation Errors
- Pial Surface
- Talairach Errors

FreeSurfer Directory Tree

Each data set has its own unique SubjectId (eg, bert)



Actual Workflow

1. recon-all --autorecon1 (Stages 1-5)
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22)
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29)


Note: all stages can be run individually

Troubleshooting: Skull Strip



brainmask.mgz

Troubleshooting: Intensity Normalization



Intensity Normalization Failure.

All WM in T1 volume (T1.mgz) should be 110.

Can fix by adding “Control Points”. Beware partial voluming!

Actual Workflow

1. recon-all --autorecon1 (Stages 1-5)
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22)
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29)



Note: all stages can be run individually

Actual Workflow

1. recon-all --autorecon1 (Stages 1-5)
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22)
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29)


Note: all stages can be run individually

Troubleshooting: Segmentation Error



Eye Socket classified as WM.
Skull Strip Failure.

Troubleshooting: Segmentation Error




“Hypo-Intensities”
White Matter Lesions

Actual Workflow

1. recon-all --autorecon1 (Stages 1-5)
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22)
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29)

Note: all stages can be run individually

Troubleshooting: Pial Surface Error



White/Gray OK, but
Pial Inaccurate

Actual Workflow

1. recon-all --autorecon1 (Stages 1-5)
2. Check talairach transform, skull strip, normalization
3. recon-all --autorecon2 (Stages 6-22)
4. Check surfaces
 1. Add control points: recon-all --autorecon2-cp (Stages 10-22)
 2. Edit wm.mgz: recon-all --autorecon2-wm (Stages 13-22)
 3. Edit brain.mgz: recon-all --autorecon2-pial (Stage 19-22)
5. recon-all --autorecon3 (Stages 23-29)

Note: all stages can be run individually

The End!

Feedback: astevens@nmr.mgh.harvard.edu